



▶ **Manual Measuring Stage MT 200 x 100 MR** for measuring in reflected or transmitted light

Measuring range:	200 x 100 mm
Accuracy:	Absolut accuracy in $\mu\text{m}$ : measuring length / 40 + 2,5
Orthogonality:	< 20 s, typ. 15 s
Resolution:	0,1 $\mu\text{m}$
Material:	Aluminium
Weight:	7.000 g
Load capacity:	5.000 g
Surface:	anodic coating, black lacquered
Certification:	DIN EN ISO 9001 : 2000
Packaging size / weight:	66 x 56 x 36 cm / 10 kg

**Order Information:**

Part No.: 00-31-104-0000

**Features:**

- ▶ for measuring in reflected or transmitted light
- ▶ for material applications
- ▶ equipped with MR measuring system
- ▶ incl. Digital Readout unit SCDplus
- ▶ single-hand control
- ▶ quick adjustment in XY
- ▶ adjustable torque

**Accessories:**

Glass plate (incl. in delivery)  
Part No.: 00-31-101-0810

For further accessories please refer to:  
[www.marzhauser.com](http://www.marzhauser.com)