



▶ **MT 50 x 50 MR** for measuring in reflected or transmitted light, with MR measuring system

	Specifications	Order Information:
Travel range:	50 x 50 mm (2" x 2")	Part No.: 00-30-119-0000
Accuracy:	absolut accuracy in μm : measuring length in mm / 40 + 2,5	Accessories:
Resolution:	0,1 μm	Metal plate
Weight:	1.400 g	Part No.: 00-30-119-0000
Load capacity::	5.000 g	Features:
		▶ for reflected or transmitted light applications
		▶ equipped with digital display SCDplus
		▶ equipped with MR Measuring System