



▶ **Manual Microscope Stage EK 100 x 100** for Upright Microscopes

<b>Travel range:</b>	max. 102 x 102 mm
<b>Drive:</b>	coaxial drive / rack and pinion
<b>Material:</b>	aluminium
<b>Surface:</b>	anodic coating, black laquered
<b>Weight:</b>	1,600 g
<b>Load capacity:</b>	5,000 g

**Order Information**

Order No.: 00-22-401-0000

**Accessory**

Glass plate (included in delivery)  
Order No.: 00-24-401-0810

For further accessories please refer to:  
[www.marzhauser.com](http://www.marzhauser.com)

**Dimensions:**

