



▶ **Manual Microscope Stage EK 100 x 80** for Upright Microscopes

Travel range:	max. 102 x 82 mm
Drive:	coaxial drive / rack and pinion
Material:	aluminium
Surface:	anodic coating, black laquered
Weight:	1,600 g
Load capacity:	5,000 g

Order Information

Order No.: 00-22-420-0000

Accessory

Glass plate (included in delivery)
Order No.: 00-24-427-0810

For further accessories please refer to:
www.marzhauser.com

Dimensions:

