



▶ **Manual Microscope Stage EK 200 x 200** for Upright Microscopes

Travel range:	max. 205 x 205 mm (6" x 6")
Drive:	coaxial drive / rack and pinion with patented coarse adjustment
Material:	aluminium
Surface:	anodic coating, black laquered
Weight:	6,500 g
Load capacity:	5,000 g

Order Information

Order No.: 00-23-801-0000

Product Features

- ▶ equipped with patented coarse adjustment
- ▶ fine adjustment via coaxial drive

Accessory

Metal plate (included in delivery)
 Order No.: 00-24-801-0011

For further accessories please refer to:

www.marzhauser.com

Dimensions:

