



▶ **Manual Microscope Stage EK 300 x 300** for Upright Microscopes

Travel range:	max. 305 x 305 mm (12" x 12")
Drive:	coaxial drive / rack and pinion with patented coarse adjustment
Material:	aluminium
Surface:	anodic coating, black laquered
Weight:	18,200 g
Load capacity:	5,000 g

Order Information

Order No.: 00-23-012-0000

Product Features

- ▶ equipped with patented coarse adjustment
- ▶ fine adjustment via coaxial drive
- ▶ for reflected light applications

Accessory

For accessories please refer to:
www.marzhauser.com

Dimensions:

