



Scanning Stage SCAN 100 x 80 for upright microscopes

Specifications

Travel range: max. 100 x 80 mm

Repeatability: < 1 µm*

measured according to VDI / DGQ 3441

Accuracy: +/- 3 µm*

measured according to VDI / DGQ 3441

Orthogonality: < 40 sec., typically < 25 sec.

Resolution: 0,05 µm* (smallest step size)

Max. travel speed: 60 mm/s*, with 1 mm lead screw pitch

120 mm/s*, with 2 mm lead screw pitch 240 mm/s*, with 4 mm lead screw pitch

Motor: 2-phase stepper motor

Material: aluminium

Limit switches: adjustable limit switches

Weight: 4.300 g

Surface: anodic coating, black
Certification: DIN EN ISO 9001 : 2000
Packaging size / weight: 58 x 48 x 38 cm / 8 kg

Order Information

Part No.: 00-24-427-0000, with 1 mm lead screw pitch Part-No.: 00-24-428-0000, with 2 mm lead screw pitch Part-No.: 00-24-429-0000, with 4 mm lead screw pitch

Features

- adjustable limit switches for individual travel range limitation
- exchangeable stage inserts for different kinds of applications
- motor cable connection front right
- optional measuring system for closed-loop operation available

Accessories See accessories list!

* Specifications valid with Märzhäuser stage controller only

Specifications subject to change without notice.

Märzhäuser Wetzlar GmbH & Co. KG In der Murch 15

D-35579 Wetzlar (Germany)

Tel.: +49 6441 9116-0 Fax: +49 6441 9116-40 eMail: info@marzhauser.com

www.marzhauser.com