



▶ **Motorized Focus Drive MFD-2** for supplementary motorization of the Z-axis

Resolution:	typical 0,002 μm Depends on the speed reduction gear ratio of the microscope's fine focus, as well as the type of controller used!
Motor:	2-phase stepper motor
Max. revolutions:	60 rev. / s
Material:	aluminium
Surface:	anodic coating, black laquered
Weight:	ca. 710 g
Dimensions:	Design and dimensions vary in accordance with the microscope it is used on.

Order Information

MFD-2 for Nikon Eclipse MA 200
Order No.: 45-54-601-0000

Accessory

The MFD-2 can be operated with any Märzhäuser stepper motor controller in conjunction with operating device Ergodrive or Joystick.

For further accessories please refer to:
www.marzhauser.com